

# APL103(Experimental Methods)

## LAB REPORT

### Exp no:- 7 (Cantilever Vibration)

1. Aditya Agrawal (2021AM10198)
2. Swapnil Kashyap (2021AM10782)

#### Objective:



1. To study damped vibrations of a cantilever beam.
2. To determine the damping ratio and un-damped natural frequency to the above system.

#### Apparatus:

- i) A cantilever beam on which a resistance-type strain gauge has been mounted near the clamped end.
- ii) A signal conditioner/strain meter and data acquisition system.
- iii) Computer for data storage and analysis.

#### Theory:

When a metal cantilever, as shown in the figure below, is given an initial displacement and released, it exhibits damped oscillations. The damping is usually very small, so the oscillations persist for a long time and can be easily recorded. In the set-up two gauges have been used, both on the top surface of the cantilever, and they are to be connected on opposite sides of the Wheatstone bridge for measuring the bending strain.



## Procedure:

- i) Start up the software and set the sampling frequency to 10 kHz.
- ii) Begin acquisition and set the cantilever in vibration by giving it an initial displacement and releasing it.
- iii) Observe the damped oscillations on the screen and stop acquiring data after the amplitude falls to less than one-tenth of the initial value.
- iv) Select the required range of data and store it in a file for analysis.
- v) Repeat the above steps four more times --- i.e., take 5 data sets in all.

## Analysis:

Using a spreadsheet such as Excel or your own software, determine the following:

- i. The time taken for 10 oscillations. Using this, determine the time period,  $T$ , and the damped frequency:

$\omega_d = 2\pi/T$  which is related to the undamped natural frequency,  $\omega_n$ , by  $\omega_d = \omega_n \sqrt{1 - \xi^2}$   
Where  $\xi$  is the damping ratio.

- ii. Pick any convenient cycle for the initial amplitude. Measure this amplitude and the amplitude after 10 cycles. Take the natural logarithm of the amplitude ratio and obtain another relation between the natural frequency and damping ratio as follows:

$$-\xi \omega_n \times 10T = \ln(A_{10} / A_0)$$

- iii. Hence, determine the un-damped natural frequency and the damping ratio.
- iv. Repeat the above steps for the other 4 data sets.

## Readings:

$T_0$	$A_0$	$T_{10}$	$A_{10}$	$T_{10} - T_0$
14.0675	0.1945	14.186	0.194	0.1185
18.9334	0.1945	19.0519	0.01939	0.1185
24.5552	0.1951	24.6746	0.1944	0.1194
30.4122	0.1944	30.5309	0.1939	0.1187
36.4445	0.1948	36.5633	0.1942	0.1188

Using the formula:  $\omega_d = 2\pi/T$

$T_{10} - T_0$	$\omega_d$
0.1185	53.02266
0.1185	53.02266
0.1194	52.62299
0.1187	52.93332
0.1188	52.88877

Now using the formula:  $\log(A_0 / A_{10}) = 2\pi\xi / \sqrt{1 - \xi^2}$

$A_0$	$A_{10}$	$\xi$
0.1945	0.194	0.000409665425
0.1945	0.1939	0.000491725182
0.1951	0.1944	0.000572059368
0.1944	0.1939	0.000409876430
0.1948	0.1942	0.000490966736

Now to calculate  $\omega_n$  we have:  $\omega_d = \omega_n \sqrt{1 - \xi^2}$

$\omega_d$	$\omega_n$
53.0226608	53.0226653
53.0226608	53.0226672
52.6229925	52.6230011
52.9333219	52.9333263
52.8887652	52.8887716

## Discussion:

1. If we take the ratio to be very small then,

$$\log ( A_0 / A_{10}) = 2\pi\xi$$

$A_0$	$A_{10}$	$\omega_n$	$\xi$
0.1945	0.194	53.0226608	0.000409665459
0.1945	0.1939	53.0226608	0.000491725241
0.1951	0.1944	52.6229925	0.000572059461
0.1944	0.1939	52.9333219	0.000409876464
0.1948	0.1942	52.8887652	0.000490966795

We can see that the error in damping ratio  $\xi$  and undamped natural frequency is negligible, i.e., the error in the average of  $\xi$  is  $5.6 \times 10^{-11}$ , and the error in the average of undamped natural frequency is  $6.06 \times 10^{-6}$ .

2. Because strain is measured along the longitudinal axis of the strip, it cannot quantify multidirectional loads. They are susceptible to overload and as a result, may be damaged. Performance is affected by humidity, temperature, hysteresis, and repeatability, and accuracy decreases with extended use.